

## **Affymetrix Probe Synthesis Guidelines AtGenExpress**

Markus Schmid and Jan Lohmann  
Department of Molecular Biology  
Max Planck Institute for Developmental Biology  
Spemannstrasse 37-39, D-72076 Tübingen, Germany

The following protocol was used for all samples contained in the developmental series of AtGenExpress (Tübingen set). The protocol is streamlined for larger sample numbers. It is imperative that you familiarize yourself with the original protocols provided by Affymetrix and by the suppliers of the individual kits before starting probe preparation.

### **RNA isolation from whole plants and seedlings**

Use of “Plant-RNeasy Mini Kit” (QIAGEN) to isolate total RNA is recommended. Make sure you use sterile and nuclease-free glassware or plastic.  
Before first use, READ THE MANUAL!!!

#### Short protocol:

- grind ~50 mg tissue under liquid N<sub>2</sub> to a fine powder
- transfer powder to an sterile Eppendorf tube
- add 450 $\mu$ l buffer RLT (containing 10- $\mu$ l/ml  $\beta$ -Mercaptoethanol) and vortex vigorously
- store sample in liquid N<sub>2</sub>; at this point you can process the remaining samples
- thaw samples
- incubate 3 min at 55°C
- apply suspension to a QIAshredder column and centrifuge for 2 min at max. speed
- transfer flow-through to a new Eppendorf tube without disturbing the pellet
- add 0.5 volumes 95% EtOH and mix by pipetting
- apply sample to RNeasy column
- centrifuge for 20 sec at 10,000g
- discard flow-through
- add 700 $\mu$ l buffer RW1 to the column and centrifuge for 20 sec at 10,000g
- transfer column to a new collection tube
- wash column with 500 $\mu$ l buffer RPE and centrifuge for 20 sec at 10,000g and discard flow-through
- repeat washing once and discard flow
- spin for an additional 2 min at max. speed to get rid of residual ethanol
- transfer RNeasy column to a 1.5ml collection tube
- elute RNA with 50 $\mu$ l RNase-free H<sub>2</sub>O; let sit on bench for 2 min before centrifugation to aid elution
- repeat elution of RNA with additional 50 $\mu$ l RNase-free H<sub>2</sub>O to a total of 100 $\mu$ l
- take OD; note concentration and 260/280 ratio
- precipitate by adding 2  $\mu$ l glycogen (Ambion, 5mg/ml), 8 $\mu$ l sodium acetate (Sigma, 3M, pH 5.2) and 250  $\mu$ l EtOH (99%) at -20°C O/N

## **RNA isolation from apices and small tissue samples**

- transfer tissue to an Eppendorf tube  
tip: To avoid losing your sample cut a small hole into a sheet of paper and place a tube in the hole. Put the whole thing over a bucket with dry ice and press the tube into the dry ice. If you spill some of the apices, just pick them up from the paper and put them into the tube
- grind tissue extensively using a eppi-pistil (work on dry ice to keep sample frozen)
- add 450 $\mu$ l buffer RLT containing  $\beta$ -mercaptoethanol
- continue grinding the tissue for 1min
- take the pistil out of the tube; make sure absolutely all tissue remains in the tube
- mix by vortexing
- store sample in liquid N<sub>2</sub>; at this point you can process the remaining samples
- thaw samples
- incubate 3 min at 55°C
- apply suspension to a QIAshredder column and centrifuge for 2 min at max. speed
- transfer flow-through to a new Eppendorf tube
- add 0.5 volumes 95% EtOH and mix by pipetting
- apply sample to RNeasy column
- centrifuge for 15 sec at 10,000g
- discard flow-through
- add 700 $\mu$ l buffer RW1 to the column and centrifuge for 15 sec at 10,000g
- transfer column to an new collection tube
- wash column with 500 $\mu$ l buffer RPE, centrifuge for 15 sec at 10,000g and discard flow-through
- repeat washing with RPE and discard flow-through
- spin for an additional 2 min at max. speed to get rid of residual ethanol
- transfer RNeasy column to a 1.5ml collection tube
- elute RNA with 50 $\mu$ l RNase-free water; let sit on bench for 2 min before centrifugation to aid elution
- repeat elution of RNA with additional 50 $\mu$ l RNase-free H<sub>2</sub>O to a total of 100 $\mu$ l
- take OD; note concentration and 260/280 ratio
- precipitate by adding 2  $\mu$ l glycogen (Ambion, 5mg/ml), 8 $\mu$ l sodium acetate (Sigma, 3M, pH 5.2) and 250  $\mu$ l EtOH (99%) at -20°C O/N

## **Resuspension of RNA**

- spin 30 min 14k rpm at 4°C to collect RNA precipitate
- carefully remove supernatant
- wash with ice-cold EtOH (70%)
- spin 15 min 14k rpm at 4°C
- carefully remove supernatant
- wash with ice-cold EtOH (70%)
- spin 15 min 14k rpm at 4°C
- carefully remove supernatant
- air dry pellet carefully (do not overdry!)
- resuspend RNA in RNase free water to obtain a concentration of at least 0.5  $\mu\text{g}/\mu\text{l}$  (maximum volume for cDNA synthesis is 10  $\mu\text{l}$ )
- use 1  $\mu\text{l}$  to take OD and ratio.
- RNA can be stored at  $-80^\circ\text{C}$  at this point

## **ds-cDNA synthesis**

Use “Superscript double stranded cDNA synthesis kit” (Invitrogen, # 11917-010) for cDNA synthesis. The kit provides all reagents needed for the synthesis of ds-cDNA. RNA spike controls are sold by Affymetrix (Poly A RNA control kit; # 900433)

### 1<sup>st</sup> strand cDNA synthesis:

Use 5  $\mu\text{g}$  total RNA in 10  $\mu\text{l}$ .

- incubate at RNA 70°C for 10 min
- chill rapidly on ice
- prepare mastermix from the following reagents for each sample:

1.1  $\mu\text{l}$  T7-dT24-oligo (HPLC purified; 100pmol/ $\mu\text{l}$ )  
4.4  $\mu\text{l}$  5x 1<sup>st</sup> strand buffer  
2.2  $\mu\text{l}$  0.1mM DTT  
1.1  $\mu\text{l}$  10mM dNTP mix  
2.2  $\mu\text{l}$  SSII reverse transcriptase  
0,22  $\mu\text{l}$  spike controls (prepare 1:20 stock, then dilute 1:50)

- add 10  $\mu\text{l}$  mastermix to each RNA sample (final volume 20 $\mu\text{l}$ ), mix gently; if necessary spin down briefly
- incubate at 45°C for 1h

### 2<sup>nd</sup> strand cDNA synthesis:

- prepare mastermix from the following reagents for each sample
  - 100  $\mu\text{l}$  H<sub>2</sub>O (RNase-free)
  - 33  $\mu\text{l}$  5x 2<sup>nd</sup> strand buffer
  - 3.3  $\mu\text{l}$  10mM dNTP mix
  - 1.1  $\mu\text{l}$  *E. coli* DNA Ligase (10U/ $\mu\text{l}$ )
  - 4.4  $\mu\text{l}$  *E. coli* DNA Polymerase I (10U// $\mu\text{l}$ )
  - 1.1  $\mu\text{l}$  *E. coli* RNaseH (2U/ $\mu\text{l}$ )
- add 130  $\mu\text{l}$  mastermix to the 1<sup>st</sup> strand cDNA (20 $\mu\text{l}$ ) on ice
- mix gently; if necessary spin down briefly
- incubate at 16°C for 2h (don't let temperature rise above 16°C)

## **ds-cDNA clean-up**

use the “GeneChip Sample Cleanup Module” from Affymetrix (# 900371), READ the manual!

- add 600  $\mu\text{l}$  cDNA binding buffer to the 150  $\mu\text{l}$  2<sup>nd</sup> strand reaction
- vortex to mix (color of mix must be yellow!)
- apply the sample mix to a cDNA Cleanup Spin Column sitting in a 2 ml collection tube
- spin at > 8000g for 1 min
- discard flow through
- add 750  $\mu\text{l}$  cDNA wash buffer (with EtOH !)
- spin at > 8000g for 1 min
- discard flow through and transfer spin column into new 2ml collection tube
- open cap and spin full speed for 2 min to remove residual EtOH
- transfer spin column into new 1.5 ml collection tube
- pipet 14  $\mu\text{l}$  cDNA Elution Buffer directly onto the membrane and incubate for 1 min
- spin full speed for 1 min to elute (average elution volume is 12  $\mu\text{l}$ )
- ds-cDNA can be stored at  $-80^{\circ}\text{C}$  at this point

## **in vitro Transcription**

The “BioArrayHighYield Transcript Labeling Kit” (Enzo; also distributed by Affymetrix) is the system of choice to synthesize biotin-labeled cRNA for use with Affymetrix GeneChips.

### Short protocol:

- thaw ds-cDNA on ice
- prepare mastermix from the following reagents for each sample
  - 11  $\mu\text{l}$  H<sub>2</sub>O nuclease-free
  - 4.4  $\mu\text{l}$  10x IVT-buffer
  - 4.4  $\mu\text{l}$  10x Biotin-labeled Ribonucleotide mixture
  - 4.4  $\mu\text{l}$  10x DTT
  - 4.4  $\mu\text{l}$  10x RNase Inhibitor
  - 2.2  $\mu\text{l}$  20x T7 RNA-Polymerase
- add 28  $\mu\text{l}$  mastermix to 12  $\mu\text{l}$  ds-cDNA (final volume 40  $\mu\text{l}$ ) mix gently; if necessary spin down briefly
- incubate at  $37^{\circ}\text{C}$  for 5h

## **cRNA clean-up**

use the “ GeneChip Sample Cleanup Module” from Affymetrix (# 900371), READ the manual!

### Short protocol:

- add 60  $\mu$ l water (RNase free) to the IVT reaction and vortex
- add 350  $\mu$ l IVT cRNA Binding Buffer and vortex
- add 250  $\mu$ l EtOH (99%) to the lysate and mix well by pipetting – do not centrifuge!
- apply sample mix (700 $\mu$ l) to IVT cRNA Cleanup Spin Column sitting in a 2ml collection tube
- spin at > 8000g for 1 min
- discard flow through
- add 500  $\mu$ l IVT cRNA Wash Buffer (with EtOH)
- spin at > 8000g for 15 sec
- discard flow through
- add 500  $\mu$ l 80% EtOH
- spin at > 8000g for 15 sec
- discard flow through and transfer spin column into new 2ml collection tube
- open cap and spin full speed for 2 min to remove residual EtOH
- transfer spin column into new 1.5 ml collection tube
- pipet 11  $\mu$ l RNase free water directly onto the membrane and incubate for 1 min
- spin full speed for 1 min to elute
- pipet 10  $\mu$ l RNase free water directly onto the membrane and incubate for 1 min
- spin full speed for 1 min to elute
- use 1  $\mu$ l (1:100 dilution in water works fine) to take OD for cRNA concentration and 260/280 ratio
- cRNA can be stored at  $-80^{\circ}\text{C}$
- calculate adjusted cRNA yield/concentration:

$$\text{adjusted cRNA yield} = \text{RNA}_{\text{m}} - (\text{total RNA} * y)$$

RNA<sub>m</sub> = amount of cRNA measured ( $\mu$ g)

total RNA = amount of total RNA used for cDNA synthesis ( $\mu$ g); usually 5 $\mu$ g

y = fraction of cDNA reaction used for in vitro transcription (usually 1)

## **cRNA fragmentation**

- mix:
  - 15  $\mu\text{g}$  cRNA (use adjusted cRNA concentration!) in a total volume of 16  $\mu\text{l}$
  - 4  $\mu\text{l}$  5x fragmentation buffer (provided in Cleanup Module)
- incubate at 94°C for 35 min
- chill rapidly on ice
- store 18  $\mu\text{l}$  fragmented cRNA at -80°C, use the remaining 2  $\mu\text{l}$  for gel analysis

## **Gel analysis**

It is a good idea to check the quality of the synthesized cRNA and the fragmentation on an agarose gel before wasting a GeneChip on a bad probe. It is not necessary to run a denaturing RNA gel, but samples should be prepared in a denaturing sample buffer\*.

### Short protocol:

- add 5  $\mu\text{l}$  sample buffer to the following samples:
  - 2  $\mu\text{l}$  cRNA
  - 2  $\mu\text{l}$  fragmented cRNA
- incubate at 65°C for 10 min
- chill on ice
- run the samples side by side on a 2% agarose gel prepared in 1x TAE, use 1kb marker (GibcoBRL) as a standard.

\* mix:

0.1  $\mu\text{l}$  EtBr (1mg/ml)

1  $\mu\text{l}$  10x MOPS

5  $\mu\text{l}$  formamide

2  $\mu\text{l}$  formaldehyde (37%)

3  $\mu\text{l}$  RNA dye (50% glycerol containing 12.5mg/ml bromphenol blue and xylene cyanole, each)

## **Preparation of Hybridization cocktail**

- prepare mastermix from the following reagents for each sample

96 $\mu\text{l}$	RNase free water
137.5 $\mu\text{l}$	2x hybridization buffer (prepared according to Affymetrix)
2.75 $\mu\text{l}$	Herring sperm DNA (10 mg/ml, Promega)
2.75 $\mu\text{l}$	acetylated BSA (50mg/ml, Invitrogen)
2.75 $\mu\text{l}$	control oligo B2 (Affymetrix)
13.7 $\mu\text{l}$	20x eukaryotic hybridization controls (Affymetrix heated to 65°C for 5 min)
- thaw the remaining 18  $\mu\text{l}$  fragmented cRNA on ice
- add 232  $\mu\text{l}$  mastermix to each fragmented cRNA for a total volume of 250  $\mu\text{l}$
- hybridization mix can be stored at  $-80^{\circ}\text{C}$

## **Hybridization**

- adjust GeneChip to room temperature
- prehybridize with 300  $\mu\text{l}$  1x hybridization buffer at 45 °C for at least 10 min
- heat hybridization mix to 99°C for 5 min followed by 5 min at 45°C
- remove prehybridization buffer from GeneChip cartridge
- fill GeneChip cartridge with 200  $\mu\text{l}$  hybridization mix
- hybridize over night at 45°C at 60 rpm

## **Washing, Staining and Scanning**

- remove hybridization mix from GeneChip cartridge and store at  $-80^{\circ}\text{C}$  for later use
- fill GeneChip cartridge completely with non-stringent wash buffer (prepared according to Affymetrix)
- wash and stain following Affymetrix protocols or store at 4°C for up to 10 h